

Aspire Rehab RX2 Tilt in Space Wheelchair



The Rehab RX2 Tilt-in-Space Wheelchair is designed to support individuals with complex seating, positioning, or pressure care needs. Its highly adjustable features and compatibility with pressure care cushions enable clinicians to tailor the chair to each user's unique requirements, facilitating effective posture management, enhanced comfort, and more efficient daily care routines.

- **Tilt-in-space up to 35°:** Allows postural variation throughout the day, facilitating pressure redistribution and enhancing comfort during prolonged seated periods
- **Lightweight aluminium frame (from 27.7 kg):** Minimises carer effort during pushing or manual handling, while supporting safer daily mobility and transport requirements
- **Customisable seating dimensions:** Adjustable seat width, depth, backrest height, armrest position, and headrest placement allow for tailored postural support to meet individual user needs
- **Contoured cushions with multi-stretch covers:** Provide stable support and help reduce pressure build-up, promoting skin protection and overall user comfort
- **Elevating and swing-away legrests:** Assist with leg positioning to improve blood circulation, support foot and ankle posture, and enable easier side or front transfers
- **Removable rear wheels and height-adjustable castors:** Enhance manoeuvrability, simplify storage in tight spaces, and support vehicle loading
- **Brake options for user and carer:** Dual systems, including push-to-lock and drum brakes, provide users with independence and carers with added safety control
- **Transport ready:** Crash-tested to AS/NZS 3691.19 standards for use in vehicle transport with appropriate restraint systems

✓ Contract approved

Compliance

AS/NZS 3695.1:2011

- **Anti-tippers included as standard:** Increase rear stability when tilting or navigating slopes and ramps

Code	Display Name	Total Width	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
MWS449780	Aspire Rehab RX2 Tilt in Space Wheelchair - 510 mm - 560 mm	670 - 740 mm	150kg	32.3 kg	450 - 580 mm	510 - 560 mm	430 - 505 mm	430 - 505 mm	410 - 730 mm
MWS449770	Aspire Rehab RX2 Tilt in Space Wheelchair - 460 mm - 510 mm	620 - 690 mm	150kg	30 kg	450 - 580 mm	460 - 510 mm	430 - 505 mm	430 - 505 mm	410 - 730 mm
MWS449760	Aspire Rehab RX2 Tilt in Space Wheelchair - 410 mm - 460 mm	570 - 640 mm	150kg	27.7 kg	450 - 580 mm	410 - 460 mm	430 - 505 mm	430 - 505 mm	410 - 730 mm