



GABRIEL X FIVE ACTION LIFT RECLINE CHAIR

The Aspire Gabriel X Five Action Lift Recline Chair is designed to support comfort, mobility, and ease of use. With powered adjustments for tilt, headrest, lumbar, backrest, and leg elevation, it helps users find and maintain a supported sitting position. Zoned massage and seat heating functions assist with relaxation and comfort.



Ebony
Large



Dark Grey
Small

2 YEAR
WARRANTY



Small



Large



Zero Gravity
Tilt-in-space



Back
Recline



Leg
Raise



Power
Lumbar



Power
Headrest



Seat
Heating



3 Function
Massage



Zero Gravity Tilt-in-space Positioning To create a sense of weightlessness to improve relaxation



Legrests Extend To provide a longer support under the lower legs, beneficial for taller people



USB A and USB C Charging: Located conveniently in the armrest so users can charge devices while browsing

Five Action: Back recline, leg raise, lumbar support, independent tilt, and headrest ideal for watching TV when reclined

3 Levels Backrest Massage: Rest with ease with a massage level that suits you

Electric Headrest Adjustment: Offers customisable neck positioning at the press of a button



Large Button Handset: Offers high visibility back-lit buttons that are easy to use

Adjustable Seat Heating: Offers comfort for users during cold temperatures

All-day Battery Backup: A lithium-ion battery ensures continued operation during power outages, supporting safety and independence

Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width
CHP227592	Dark Grey	Petite	475mm	450mm	490mm	710mm	1100mm	850mm
CHP227597	Ebony	Large	500mm	520mm	530mm	780mm	1170mm	900mm