



Genesis HD Adjustable Seat Walker

The Aspire Genesis HD is extra strong to support users weighing up to 204kg, yet it is still lightweight and easy to lift and manoeuvre due to the strengthened aluminium elliptical tube frame. It is highly adaptable and comfortable with adjustable seat and handle heights, padded seat and generously wide padded backrest. Rugged solid 8” wheels are ideal for both indoor and outdoor terrain.

Adjustable Seat Height: Suitable for variable user sizes, optimising seat height for the user assists with sitting and standing easily and safely

Adjustable Handle Height and Soft-touch Brakes: Correct handle height promotes walking posture while soft brakes are nice and easy to operate

Wide Padded Seat and High Weight Capacity: Ensures comfort and stability for larger users, with the added benefit of a discreet storage pocket beneath the seat

Large Padded Backrest: Provides generous and comfortable lumbar support when seated

Oversized Under-seat Bag: Ample storage space for personal items, accessible by lifting the seat

Easy Folding Mechanism and Lightweight: Lift to fold for compact storage and easy handling at only 11.3kg total weight

Solid 8” Wheels with Inline Rear Brakes: Tough treaded wheels for a smooth ride in variable terrain, while inline rear brakes maximise manoeuvring space for the user



Marble Blue

 **1 YEAR LIMITED WARRANTY**

 **SWL 204kg**

| Code | Overall Height mm | Overall Length mm | Overall Width mm | Seat Width mm | Seat Height mm | Width Between Arms mm | Castor Size mm | Product Weight kg | SWL Under Seat Bag kg |
|-----------|----------------------|----------------------|---------------------|------------------|-------------------|--------------------------|-------------------|----------------------|--------------------------|
| WAF750336 | 875 - 1050 | 700 - 800 | 660 | 450 | 500-590 | 520 | 200 | 11.3 | 5 |

Note: Range of motion dimensions are accurate to a +/-10mm tolerance