



PACE FOREARM WALKER

The Aspire Pace Forearm Walker is designed to support balance and stability while helping the user to walk with an upright gait. This unique aid is ideal for individuals with limited hand and wrist strength as lower arm supports allow them to bear upper body weight comfortably and confidently.



Ergonomic Handles with Brakes

Depth and angle adjustable to ensure optimum comfort and accessibility to brakes and handles for control and support.



Height-Adjustable Design

The soft padded lower arm support platform (including brakes & handles) has lever activated height adjustment to achieve the optimum height for the individual at any time, and to adjust as needs change. Includes gas assisted raise and lockable height setting once adjusted.



Intuitive Braking System

Easily activated brakes provide intuitive control of rear castors to assist support and balance while walking.



Swivel Castors

All castors are total locking to ensure the forearm walker can remain safely stationary and stable when required.



1 YEAR WARRANTY

**SWL
150kg
23.6st**

	Overall Length	Overall Height	Overall Width	Width Between Forearm Supports	Wheel Size	Product Weight
Pace Forearm Walker						
	800mm	915 – 1325mm	640mm	660 – 790mm	100mm	18kg
WAF699220	31.5"	36 – 52"	25"	26 – 31"	4"	39.7lbs

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk