




SHOWER STOOL WITH ARMS



BTS119000



BTS119025

- Clip on seat or padded back attachment for easy assembly and cleaning
- Easy push button height adjustable legs
- All frames individually treated for corrosion resistance
- Seat with anterior tilt to facilitate standing from sitting and provides greater comfort for stiff or painful hips
- Raised armrests for assistance in transferring on and off the stool

- Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
- Additional height increments on one leg to accommodate uneven surfaces

Tested to AS/NZS3973:2009

DVA Contracted

Registered in the ARTG, AUST R ARTG 230166/230167

5 YEARS WARRANTY
Corrosion

1 YEAR WARRANTY
Structure



Shower Stool With Arms

	Overall Height mm	Overall Width mm	Overall Depth mm	Seat Height mm	Seat Width mm	Seat Depth mm	Between Arms mm	Product Weight kg
Steel BTS119000	670 - 855	545	400	505 - 655	575 - 775	310	495	3.8
Aluminium BTS119005	670 - 855	545	400	505 - 655	575 - 775	310	495	3.1
With Padded Seat BTS119020	670 - 855	545	400	505 - 655	575 - 775	310	495	4.2
With Padded Seat & Back BTS119025	805 - 955	545	400	505 - 655	575 - 775	310	495	6.8