



VIKING FOREARM WALKER

The Viking Forearm Walker is designed for hospital and clinical use, providing stable, comfortable, and fully adjustable support for upright walking. Its strong aluminium frame, padded arm supports, and easy-to-use hydraulic height adjustment make it ideal for rehabilitation and assisted mobility in care settings.



- **Hydraulic Height Adjustment:**
Lever-operated system makes it easy to raise or lower the arm supports smoothly and safely, helping users of different heights
- **Adjustable Padded Arm Supports:**
Soft armrests offer comfortable forearm support, for users with limited grip strength. The width can be adjusted to suit individual body shapes and provide added stability
- **Width-adjustable Base:** Improves manoeuvrability in tight spaces and around clinical equipment
- **Angle-adjustable Soft Handgrips:**
Supports a natural wrist position, making it easier and more comfortable for users to maintain control
- **Oxygen Tank Holder:** Included as standard, which can be fitted or removed to suit the user's needs
- **Optional Accessories:** Available to enhance functionality

1 YEAR WARRANTY



	Overall Length	Overall Height	Overall Width	Width Between Forearm Supports	Wheel Size	Product Weight	Bag SWL
Viking Forearm Walker							
WAF699130	850mm	900 - 1325mm	730 - 890mm	520 - 750mm	100mm	18.3kg	5kg