



Vogue Carbon Seat Walker

Carbon Fibre Frame Extremely lightweight to maximise portability

One of the lightest seat walkers on the market. Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super-lightweight carbon fibre frame for a superior walking experience.

Less Than
5kg
11lbs




Deep Seat and Wide Backrest
For improved comfort and support




Compact Turning Circle
For improved mobility in tight spaces




1 YEAR WARRANTY



Sapphire
Short



Black
Medium / HD



Emerald Green
Tall





Ergonomic Handles

These offer a comfortable grip during extended use. The flat-top design provides a stable resting place for arms while seated



Soft-touch Brakes

Lock the rear wheels to enhance safety and stability during use



Adjustable Handles

Allows for personalised positioning, helping to maintain proper posture



Simple Folding Mechanism

Locks for easy transportation and storage



Removable Storage

Zip bag secures up to 5kg of essentials



Accessories

With a cane holder included as standard, choose from our range of optional accessories

	Overall Length	Overall Width	Overall Height	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight*
	mm	mm	mm	mm	mm	mm	mm	kg	kg
Sapphire - Short WAF705340	655	600	715 - 830	460	500	230	203	4.7	5.9
Black - Medium WAF705350	655	615	725 - 940	460	550	230	203	4.9	6.2
Emerald Green - Tall WAF705370	655	600	780 - 1030	460	620	230	203	5.2	6.3
Black - HD WAF705375	655	690	780 - 1030	550	620	230	203	5.9	6.6

AS/NZ Approved | Vogue is registered in the ARTG, AUST R ARTG 271255 | *Includes accessories, bag, and backstrap.