



ZULU SCOOTER

MAXI

Powerful, Durable, and Built for the Long Journey

The Aspire Zulu Maxi mobility scooter is designed for larger users or longer distances. With a travel range of up to 50km per charge, robust suspension, and enhanced safety features, this heavy duty scooter is ideal for users seeking comfort, confidence, and independence in outdoor settings.





Reliable Performance

High-torque 1000W motor also features 1700W peak power to turbo up steeper inclines.



Digital Dashboard

Clearly shows battery life, speed, and other key information, giving you confidence while you travel.



Comfortable Design for Larger Users

The adjustable, large padded seat provides comfort and assists easy transfers. Stable platform supports a 226kg safe working load.



Robust Outdoor Performance

Front and rear wheel suspension and electro-mechanical brakes ensure control and stability on uneven terrain.



Long-distance Capability

110Ah batteries support up to 50km of travel per charge, reducing the need for frequent charging.



Smart Safety

Emergency stop button and automatic speed reduction while turning keeps you safe and secure.

Overall L x W x H	1600 x 720 x 1280mm
Seat Size D x W	510 x 560 - 610mm
Seat to Deck Height	460 - 510mm
Seat to Ground Height	670 - 720mm
Armrest Height	180 - 240mm
Min Turning Radius	1860mm
Ground Clearance	110mm
Max Climbing Angle	9°
Max Driving Speed	10km/h
Driving Range	50km
Tyres	380mm Pneumatic
Battery	2 x 12V 110Ah (Lead Acid)
Product Weight (excl. battery)	112.7kg
Battery Weight	65.3kg
Power (Nominal Peak)	1000W 1700W
Motor	4 Pole
Controller	RHINO 2 120A
Brake System	Electro-Mechanical
Free Wheel Mode	Yes
Charger Size	24V 8A

*Scooter Frame 7 Years | Motors/Axles/Brakes 5 Years | Electronics 1 Year | Battery 6 Months



SCT671320
Midnight Blue



020426