

ComfiMotion® Collection

Luxurious adjustable bed collection for an enhanced sleeping experience




 aspire

Activ Care Bed

Ultra Low Adjustable

Function, Style & Support



Supporting the modern, independent lifestyle.

The Activ Care bed combines the best of in-home care, functionality, positioning, and low height safety with luxurious and contemporary style for enhanced comfort, care, and wellbeing.

- **Ultra-low Height:** Minimum platform height of 195mm reduces falls risk, improving user safety.
- **Customisable Safe Entry/Exit Height:** Programmable single button press ensures perfect bed height to maximise independence.
- **Underbed Clearance:** Provides clear access for patient lifters and overbed tables.
- **Battery Backup:** Provides confidence and maintains independence during power outages.
- **Night Assist:** Luminous underbed halo for safe entry and exit.
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder, and back relief.
- **Pull-apart Design:** Enables more accessible delivery into homes, aged care and other health facilities.

Functions

- Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks.
- Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring, and sleep apnoea.
- Auto-regression in the Backrest:** Stops the user sliding down the bed during movement while improving comfort.
- Clinical Positioning:** One touch raises back and legs to alleviate cardiac exertion, improving resting comfort.

- Leg Raise:** Promotes circulation, relieving tired, swollen, fatigued, or painful limbs.
- Memory Function:** Allows the user to set and return to their preferred position at any time with one click of a button.
- Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- Seated Position:** Assists with bed entry and exit and upright dining.
- Trendelenburg (and Reverse):** Can assist care givers when repositioning a user*.
- Lockable Castors:** Eight heavy duty dual castors.



Companion Sleeping
See Configurations section for more companion options.

Sizes Length x Width x Height Range in mm


Long Single¹ 2030 x 920 x 195 - 700 BEB046745	King Single¹ 2030 x 1070 x 195 - 700 BEB046755	Long Double 2030 x 1370 x 195 - 700 BEB046765	Queen 2030 x 1520 x 195 - 700 BEB046775
---	--	--	--

- 10 YEAR WARRANTY
Frame
- 2 YEAR WARRANTY
Electronics
- 1 YEAR WARRANTY
Handset and Cables



Care Bed

Variable Height Adjustable



For comfort you deserve.

The Care adjustable bed blends seamlessly into the home without compromising on the therapeutic needs of the user or safety of partners or carers. Select from multiple position options designed to maximise independent mobility and improve sleep quality.



Functions

- Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks.
- Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring, and sleep apnoea.
- Clinical Positioning:** Raises back and legs to alleviate cardiac exertion, improving resting comfort.
- Lockable Castors:** Four heavy-duty castors.
- Leg Raise:** Promotes circulation, relieving tired, swollen, fatigued, or painful limbs.
- Seated Position:** Assists with bed entry and exit and upright dining.
- Trendelenburg (and Reverse):** Can assist care givers when repositioning a user*.

Sizes Length x Width x Height Range in mm
Long Single¹ 2030 x 920 x 280 - 5802 **BEB046720**
King Single¹ 2030 x 1070 x 280 - 5802 **BEB046730**

Features

- Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- Illuminated Handset:** Large button interface is intuitive and highly visible at all times.
- Enclosed Frame Design:** Ensures optimal positioning when getting into bed and safely covers the moving base platform.
- Safe-lowering Technology:** Ensures that back, leg, and Hi-Lo functions only lower under the force of gravity, reducing potential for accidents or injuries.
- Box-on-box Design:** Facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- Variable Bed Height Options** in mm
Low Height 280 - 580 with Low Profile Bed Legs
Standard 350 - 650 with Standard Castors
Accessibility 400 - 700 with Large Castors optional

10 YEAR WARRANTY
Frame

2 YEAR WARRANTY
Electronics

1 YEAR WARRANTY
Handset and Cables




¹Fixed Companion Bed Available

²With Low Profile Bed Legs.

Relax Bed

Fixed Height Adjustable



Supporting the modern, independent lifestyle.

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.



Functions

- Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring, and sleep apnoea.
- Leg Raise:** Promotes circulation, relieving tired, swollen, and/or painful limbs.
- Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button.

Sizes Length x Width x Height Range in mm
Long Single¹ 2030 x 920 x 280 - 700 **BEB046900**
King Single¹ 2030 x 1070 x 280 - 700 **BEB046910**
Long Double 2030 x 1370 x 280 - 700 **BEB046920**
Long Queen 2030 x 1520 x 280 - 700 **BEB046930**

Features

- Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- Massage:** Relaxing vibration modes for personalised neck, shoulder, and back relief.
- Press and Hold Safety Handset:** Convenient remote control handset enhances usability.
- Box-on-box Design:** Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- German Engineered Motors:** For smooth, quiet adjustment.
- Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
- USB Charge:** Connect and charge your devices conveniently.

10 YEAR WARRANTY
Frame

2 YEAR WARRANTY
Electronics

1 YEAR WARRANTY
Handset and Cables



*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques. Visit confimotion.com.au for Limited Warranty T&Cs

Mattresses & Accessories



Memory Foam Mattress (Length x Width x Height in mm)

Contoured Immersion: Layered memory foam provides support and pressure relief.

Triple-layered Foam: Top: Immersive visco-elastic layer, Middle: High resilience foam for support, and pressure redistribution. Base: Strong support for maximum durability.

Integrated Side Walls: Increase roll resistance at the bed edge.

Superior Breathability: Up to 8x more breathable than conventional memory foams.

Anti-microbial: Environmentally safe treated foam which reduces odours and kills 99.9% of bacteria.

Eco-friendly GECA Foam: Produced using cleaner chemicals and production technology.

Pressure Care Incontinence Cover: Breathable, 4-way stretch, and waterproof.*

5 YEAR WARRANTY
Foam

2 YEAR WARRANTY
Cover

Classic

Long Single
2030 x 920 x 200 **BEM469910**

King Single
2030 x 1070 x 200 **BEM469920**

Long Double
2030 x 1370 x 250 **BEM469922**

Queen
2030 x 1520 x 250 **BEM469924**

*Pressure Care

Long Single
2030 x 920 x 200 **BEM469930**

King Single
2030 x 1070 x 200 **BEM469940**

Long Double
2030 x 1370 x 250 **BEM469942**

Queen
2030 x 1520 x 250 **BEM469944**



Pressure Care



Pocket Spring Mattress (Length x Width x Height in mm)

Transitional Layered Pocket Springs: Arranged in weight activated support zones.

Integrated Side-walls: Create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.

Luxuriously Soft Pillow: Top layer provides a plush immersive feeling complimented by the responsive pocket springs.

Long Single
2030 x 920 x 240 **BEM469800**

King Single
2030 x 1070 x 240 **BEM469810**

Long Double
2030 x 1370 x 240 **BEM469820**

Queen
2030 x 1520 x 240 **BEM469830**

5 YEAR WARRANTY


Breeze Pillow Cooling Gel Memory Foam

Adaptive and supportive sleep surface
Helps regulate sleep surface temperature
Accommodates side and back sleepers
BEP050100



Plush Pillow Memory Foam and Fibrefill

Familiar sleep surface
Improved head and neck support
Accommodates side and back sleepers
BEP050110



Allure Headboard Free Standing

Stylish freestanding headboard. Subtle square padded design.


Long Single **BEA830120**

King Single **BEA830130**

Long Double **BEA830132**

Queen **BEA830134**

King **BEA830136**



Chesterfield Headboard Free Standing

Stylish freestanding headboard. Chesterfield padded design.

Long Single **BEA830100**

King Single **BEA830110**

Long Double **BEA830112**

Queen **BEA830114**

King **BEA830116**



Allure Headboard & Footboard Fixed

Stylish pivoting headboard moves with bed to support head during tilt. Footboard keeps bedding in place during profiling and can support an air mattress pump.

	Headboard	Footboard
Long Single	BEA830150	BEA830160
King Single	BEA830152	BEA830162
Long Double	BEA830154	BEA830164
Queen	BEA830156	BEA830166



Bed Rail

Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.

Care and Relax Bed **BEA830210**

Activ Care Bed **BEA830212**



Crook Handle

Provides multiple grip points to assist repositioning in bed, and safe entry and exit.

Care and Relax Bed **BEA830200**

Activ Care Bed **BEA830211**



Extendable Side Rail

Provides a telescoping bed rail support for users at risk of falls.

Activ Care Bed **BEA830205**



Split King Upgrade Kit

Ensures simultaneous Hi-Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.

Activ Care Bed LS and KS only
BEA830250



Simplified Handset

For clients with reduced cognitive ability.

Activ Care Bed **BEA830245**



Partner Handset

Allows double and queen beds to be operated from either side of the bed.

Activ Care Bed **BEA830240**



Large Accessibility Castors

Larger castors raise the bed frame to 155mm, facilitating overbed table and patient lifter access.

Care Bed only (4 x required)
MWT451080



Battery Backup Replacement

Provides peace of mind in the event of a power outage.

Care and Relax beds **BEA830220**

Activ Care beds **GSP652706**



Clinical Assistance Guide

Scan for product specifications, assessment tools, and prescription support.



AIDACARE PTY LTD

Building 3A, 1 Moorebank Avenue,
Moorebank NSW 2170
Australia

1300 133 120 | aidacare.com.au

AIDACARE LIMITED

Arcadia House, Maritime Walk,
Ocean Village, Southampton SO14 3TL
United Kingdom

01622 541235 | aidacare.co.uk

MADE IN PRC

ASBR28.102 151025