

SWEP SCRIPT FORM

OFF CONTRACT POWER MOBILITY RANGE

REFERRER DETAILS		CLIENT DETAILS	
Name:	Position:	First Name:	Last Name:
Phone:	Fax:	Aidacare Representative:	
Email:	Freight: Metro \$65+GST Regional \$80+GST		
PLEASE FORWARD FORM TO SWEP – FIRST \$4000 IS PAID BY SWEP – CO-CONTRIBUTION BELOW IS PAID BY CLIENT			



Black
black arm covers as standard

WHILL MODEL C2

Description	PWS655400 - Power Wheelchair - 450mm PWS655410 - Power Wheelchair - 500mm		
Warranty	Base & Seat Frame: 5 Years Electrical & Mechanical Components: 2 years Batteries & Wheels: 6 months		
Black	PWS655400	Seat Depth	400mm
Co-Contribution	\$0.00	Maximum Speed	8km/h
Overall Width	604mm	Backrest Height	340-460mm
Seat Height	405 – 485mm	Range	18km
Seat Width	450mm	SWL	136kg
Overall Height	985mm		



Carbon Fibre

PRIDE LITESTREAM FOLDING CARBON FIBRE

Description	Lightweight		
Warranty	Frame & Electrical Components: 2 years Battery: 1 Year		
Black	PWS654185	Maximum Speed	6km/h
Co-Contribution	\$0.00	Seat to Footrest	475mm
Total Length	930mm	Range	15km
Total Width	577mm	Material	Carbon Fibre
Seat Depth	430mm	SWL	136kg
Seat to Floor Height	475mm		



Red

PRIDE JAZZY SELECT ELITE HD

Warranty	2 Years (limited)	Total Width	616mm
Red	PWS654250	Maximum Speed	6.5km/h
Co-Contribution	\$0.00	SWL	204 kg
Total Length	920mm		



Storm
Grey



Desert
Sand



ASPIRE JUNO

Description	4 Wheel		
Warranty	Scooter Frame/Motors/Axles/Brakes: 2 Years Battery 6 Months: Electronics 1 Year		
Storm Grey	SCT671296	Total Height	1400mm
Desert Sand	SCT671295	Maximum Speed	10km/h
Co-Contribution	\$0.00	Unit Weight	147.5kg
Total Length	1490mm	Range	40km
Total Width	720mm	SWL	181kg



Silver



ASPIRE ZULU

Description	4 Wheel		
Warranty	Scooter Frame: 3 Years Motors / Axles / Brakes: 18 Months Battery: 6 Months		
Silver	SCT671310	Maximum Speed	10km/h
Co-Contribution	\$1750	Unit Weight	178kg
Total Length	1600mm	Range	50km
Total Width	720mm	SWL	226kg
Total Height	1280mm		



White



ASPIRE SUPALITE

Description	4 Wheel		
Warranty	Scooter Frame / Motors / Axles / Brakes: 2 Years Battery: 6 Months		
White	SCT670616	Maximum Speed	6.5 km/h
Co-Contribution	\$0.00	Unit Weight	34kg
Total Length	1022mm	Range	18.5km/h
Total Width	480mm	SWL	120kg
Total Height	880mm		



Red



Yellow



PRIDE PATHRIDER 130XL PURSUIT

Description	4 Wheel		
Warranty	2 Years, 6 months - Batteries		
Red	SCT670571	Total Width	635mm
Yellow	SCT670575	Maximum Speed	10km/h
Co-Contribution	\$0.00	Range	25km
Total Length	1360mm	SWL	159kg



Blue



PRIDE PATHRIDER 10 DXE

Description	4 Wheel		
Warranty	2 Years, 6 months - Batteries		
Blue	SCT670507	Maximum Speed	10km/h
Co-Contribution	\$0.00	Range	39km
Total Length	1194mm	SWL	181kg
Total Width	560mm		