

# KINERIS PRESSURE REDUCING CUSHION RANGE

Ideal for patients at very high risk of pressure injury development. It is suitable for both active and non-active users, accommodating a large range of user heights and weights.



- **Accurate Pressure Adjustment:** Manometer enables exact pressure adjustment every time. Easy for carers!
- **Single-valve:** For simple inflation and maintenance.
- **Reduced Friction and Shearing:** Cell mobility minimises friction and shearing forces.
- **Optimal Height:** Available in low and high profiles.
- **Comfortable and Breathable:** The removable cover is stretchable, and breathable, and aids body temperature regulation during prolonged sitting.
- **Stays in Place:** The inner cover features a non-slip base to keep the cushion secure.
- **Versatile Sizing:** Available in various sizes to accommodate different needs.

Code	Description	Profile	Width x Depth	Height
PTC627100	Air Cushion	Low	16.5 x 16.5in	2.75in
PTC627110	Air Cushion	Low	16.5 x 18.5in	2.75in
PTC627120	Air Cushion	Low	18.5 x 18.5in	2.75in
PTC627250	Replacement Cover	Low	16.5 x 16.5in	
PTC627255	Replacement Cover	Low	16.5 x 18.5in	
PTC627260	Replacement Cover	Low	18.5 x 18.5in	

Code	Description	Profile	Width x Depth	Height
PTC627200	Air Cushion	High	16.5 x 16.5in	4in
PTC627210	Air Cushion	High	16.5 x 18.5in	4in
PTC627220	Air Cushion	High	18.5 x 18.5in	4in
PTC627265	Replacement Cover	High	16.5 x 16.5in	
PTC627270	Replacement Cover	High	16.5 x 18.5in	
PTC627275	Replacement Cover	High	18.5 x 18.5in	