

Aspire PosturePro Plus Chair



A Practical Solution for Complex Postural and Comfort Needs. The PosturePro Plus is designed for individuals with complex or evolving needs, providing high levels of postural support and effective pressure care. It combines memory foam, air bladders, and multi-stretch fabric to deliver targeted pressure relief while maintaining user comfort.

- **Customisable Backrest - Adjustable Head, Neck, and Shoulder Sections:** Provides support for kyphotic or complex fixed postures, while accommodating contractures with a negative angle leg rest and height-adjustable armrests.
- **Integrated Laterals:** Offers lateral support to promote upright posture and prevent leaning.
- **Adjustable Pocketed Waterfall Backrest Design:** Provides customised back support through the use of rear-side pockets that allow internal padding to be added or removed. This enables the support surface to adapt to the user's specific shape and postural needs.
- [Watch the instructional adjustment video.](#)
- **Tool-Free Seating Adjustment - Seat Depth and Width Configuration:** Allows therapists and carers to quickly adjust seat dimensions without tools, ensuring individualised comfort for each user.
- **Adaptable Footrest - Footrest Depth and Angle Adjustability:** Allows customisation of footrest depth and angle to accommodate the user's positioning and weight-bearing requirements, aiding in sit-to-stand transfers.
- **Stand Air Compatibility - Safe Standing Transfers:** Integrates seamlessly with the Aspire Star Mover to ensure safe and efficient standing transfers, promoting user mobility and independence.
- **Arm Height Adjustability - Height and Depth Adjustable Armrests:** Allows adjustment of armrests to the correct height and depth, improving arm

✓ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

positioning for comfort and support.

- **Tilt-in-Space - Infinite Positioning Options:** Enables tilting and reclining at various angles, making repositioning easier and enhancing comfort.
- All features and supports of the PosturePro Plus can be adjusted over time, providing increasing support as needed for individuals with degenerative conditions
- Frame - 2 Years
- Components - 1 Year
- Fabric - 1 Year

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP231610	Aspire PosturePro Plus Chair	888 mm	1111 - 1400 mm	140 kg	450 - 540 mm	330 - 530 mm	530 - 590 mm	530 - 590 mm	777 - 900 mm