

Aspire Oxygen Bottle Holder



- 140 mm diameter stainless steel holder suits sizes Oxygen Bottle sizes C to D
- Adaptable clamp fits most shower commode brands

Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Material
BTC067200	Oxygen Bottle Carrier	Stainless Steel