

Aspire Lifecomfort Bariatric Mattress



The LIFECOMFORT® BARIATRIC is a robust support surface designed for larger users up to 300 kg. Four foam tiers with a visco-elastic memory foam upper provide immersion for pressure reduction and comfort as well as support. Side wall foam density is strengthened to assist heavier users for bed entry and exit. Fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care, Aged Care and up to Acute Care environments.

- **Durable Hospital Grade Foam** that is treated for antimicrobial and antibacterial properties
- **Breathable, 4-Way Stretch Waterproof Cover** to reduce shear, friction and maintain skin microclimate
- **Medical Grade Waterfall Cover** is fully-welded over the zipper to maximise infection control
- **Asset Management Ready** with unique batch and serial numbers printed on all covers
- **Carry Handle** - Integrated handles for safe manual handling
- **Hospital Grade** - Suitable for bleach concentrations up to 10,000 ppm
- **4-sided Zipper** - Easily removed for effective cleaning
- **4-way Stretch** - For shear mitigation and friction reduction
- Immersive surface contours and cradles around the user for personalised pressure relief and comfort
- Machine Washable Covers: Up to 95°C temperature
- **Foam Layers:**
 - IMMERSIVE SURFACE - Immersive visco-elastic
 - COMFORT LAYER - Transition layer that boosts visco layer efficiency under higher weight loads
 - SUPPORT LAYER - High resilience foam for diverse user weight and

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_043306003_0103_1_2

Power Bed –

Composite:05_423_0103_1_2

Click [here](#) for more information

Compliance

BS7177:2008

- distribution
- SAFETY & SUPPORT BASE - Strong base support for maximum durability
- INTEGRATED SIDE WALLS - Strong rigid side layers for transfer efficiency and safety
- **Applications** - Pressure Care: Medium - High
 - Hospital
 - Aged Care
 - Home Care
 - Palliative Care
 - Bariatric Care

Warranty Information

- 5 yrs Foam
- 5 yrs Cover

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM046700	Lifecomfort Bariatric Mattress - King Single	2030 mm	1050 mm	200 mm	150-300 kg