

Aspire Commode Bowl with Handle & Lid



- Aspire Commode bowl and lid with handle for transportation
- Autoclavable for sanitising
- Moulded plastic construction
- Compatible with Aspire Shower and Bedside Commodes
- Diameter: 280 mm
- Height: 160 mm

Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name
BTC081610	Aspire Commode Bowl with Handle & Lid
GSP491470	Aspire Bariatric Shower Commode - Commode Bowl