

TabTimer 7 Alarm Pill Box

TabTimer 7 Alarm Pill Box

Perfect companion for carrying daily pills with you



- Up-to 7 alarms per day
- 7 compartments for storing different medications
- Audible alarm, Vibration reminder, Light reminder, Missed pill reminder
- Easy to set and use
- When due, alarms sounds 1 minute and vibrates
- Red light flashes in compartment indicating corresponding compartment for the alarm
- Short reminder bell and lights flashing every three minutes for 15 times (45 mins)
- Press any button to stop
- 12 hour / 24 hour clock (military style)
- LCD display
- LED Backlight
- Volume level approx 60-70db
- Requires 2 x AAA size battery
- Materials ABS + PCTG
- Unit Dimensions: approx. 200 x 91 X 23mm
- Each compartment is approx. 35mm x 23mm x 15mm
(large enough for approx. 3-4 large vitamin tablets, 7-8 tablet capsules)

Code	Display Name
------	--------------

DLG275070	TabTimer - 7 Alarm Pill Box
-----------	-----------------------------

Code	Display Name
------	--------------