

JAY - Care Backrest



- Designed to accommodate a client who has a fixed kyphotic, the JAY-CARE backrest provides solid back support to accommodate and slow-down back deformity
- The contoured design and 39° angle adjustability help further accommodate the multiple needs of the kyphotic client and make the JAY-CARE suitable to a client who has changing postural needs over time
- Provides mild to moderate lateral support through the built-in lateral contour of the shell
- Optional foam wedges provide for a more customised fit
- All models include: Jay AirExchange cover and universal mounting hardware

Warranty Information

- 2 Years Backrest
- 6 Months Cover

✓ Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181003192_0105_1_2

Wheelchair Manual – Composite:

Wheelchair Powered –

Composite:05_417_0105_1_2

05_419_0105_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Safe Working Load	Seat Width	Backrest Height	Backrest Width
PTB532530	JAY Care - 16in	113 kg	15 - 16 in	16.5 in	16 in

Code	Display Name	Safe Working Load	Seat Width	Backrest Height	Backrest Width
PTB532520	JAY Care - 18in	113 kg	17 - 18 in	17.25 in	18 in
PTB532515	JAY Care - 20in	113 kg	19 - 20 in	19 in	20 in