

Elevating Legrest



- A great range of accessories to fit the Aspire Manual Wheelchair range to improve functionality and user experience
- Click [here](#) to ensure this accessory fits your specific wheelchair model

Warranty Information

- 1 Year



✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

Compliance

AS3696-1:2008

Code	Display Name
WCA848100	Elevating Legrest - Left - 300/320/350/400/410mm
WCA848110	Elevating Legrest - Right - 300/320/350/400/410mm
WCA848120	Elevating Legrest - Left - 450/460/500/510mm
WCA848130	Elevating Legrest - Right - 450/460/500/510mm
WCA848140	Elevating Legrest - Left - 570/610mm
WCA848150	Elevating Legrest - Right - 570/610mm

Code	Display Name
WCA500070	Elevating Legrest - Left - 450 to 610mm - Assist 3, Transit 3
WCA500075	Elevating Legrest - Right - 450 to 610mm - Assist 3, Transit 3
WCA500080	Elevating Legrest - Left - 400mm - Assist 3, Transit 3
WCA500085	Elevating Legrest - Right - 400mm - Assist 3, Transit 3