

Malvestio Botero Static Mattress



- Two layer foam design with high density base and castellated surface caters to bariatric patients
- Multi-stretch cover with welded seams for stringent infection control
- Strong high density side walls improve patient transfer efficiency while cut back hinge points allow the mattress to easily move with the mattress platform

Warranty Information

- Foam: 5 Years
- Cover: 2 Years



Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Compliance

BS7177:2008

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM050200	Botero Static Mattress	2000 mm	1220 mm	200 mm	400 kg