

Malvestio VIVO ICU Bed



The VIVO bed was designed with patient and care staff welfare in mind. It offers integrated features to ease care staff burden whilst creating a comfortable and effective patient support surface.

- Lateral tilting patient surface
- Alternating air mattress with full bed integration to maximise therapeutic benefit
- Patient Anti Slipping System (PAS) offsets the patient movement towards the foot end of the bed during backrest adjustment
- Radio transparent patient surface for fluoroscopic examinations
- Touch screen panel in the foot end allows for the bed and integrated air mattress functions to be controlled with ease
- Integrated bed extension to suit a range of patient heights
- Intuitive control panels within the gas assisted folding side rails
- Safety alarm system ensures the base of the bed is clear before any movements, preventing entrapment and damage
- Backrest angle of 0 to 68 degrees
- Knee bend Angle of 0 to 25 degrees
- Trendelenburg - 15 to 15 degrees

Warranty Information

- 2 Years



Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB050100	Malvestion VIVO ICU Bed	2190 mm - 2420 mm	1080 mm	480 - 820 mm	260 kg